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## 1 - Description and intended use

A partial or total shoulder prosthesis is comprised of various components, all of which are packaged individually and are designed to replace all or part of the natural shoulder joint.

A total or partial shoulder arthroplasty is intended to reduce pain and improve the mobility of the patient.

The performance capabilities of the device and its service life are closely linked to a multitude of factors, which are dependent upon the patient (physiology, occupation) and upon the surgery undertaken by the surgical team (selection, fitting and fixing of components).

## 2 - Materials

Abbreviation	Designation	Reference standard
Ti6 Al4 V-ELI	Titanium, aluminum and vanadium alloy	ISO 5832-3
UHMWPE	Ultra-high molecular weight polyethylene	ISO5834-1 & 2
T40	Titanium (unalloyed)	ISO 5832-2

## 3 - Indications

### Hemiarthroplasty with intact glenoid:

- Post-fracture necroses
- Aseptic osteonecroses
- Centered omarthroses

### Total anatomical prostheses:

Non-cuff-deficient arthropathies:

- Centered primary and secondary omarthrosis
- Inflammatory arthropathies
- Osteonecroses (ONA)
- Post-fracture malunions

### Total reverse prostheses:

- Omarthrosis/cuff-deficient inflammatory arthropathy
- Rotator cuff arthropathy
- Massive (irreparable) tear of the rotator cuff, with migration and without arthrosis
- Eccentric omarthrosis
- Post-fracture malunions
- Traumatic sequels
- Comminuted fractures in very elderly patients
- Chronic luxation of the shoulder (anterior or posterior)
- Revision of a hemiarthroplasty or a total anatomical prosthesis
- Fractures – luxations of the shoulder

## 4 - Counter-indications

The implantation of a total or partial shoulder prosthesis, in common with any surgical procedure, has attendant risks, specifically associated with anaesthesia and with the operation itself.

### General risks:

- Systemic or local infections,
- Severe vascular or neurological deficiencies
- Bone loss or poor bone quality, which is liable to compromise the stability of the implant
- Allergies to materials present in the implant
- Bone tumour in the region of the implant

### Total anatomical prostheses:

- Severe deficiency of the cuff
- Insufficient bone stock for the attachment of implants
- Permanent instability

## 5 - Warnings and precautions

Prior to any clinical use, the surgeon must be fully versed in all aspects of the operation, the implants and associated instruments, which are specially designed by MOVE-UP, together with their applications and limits. Surgeons must notify their patients of the limitations of the prosthesis, and of the associated impact upon their activities (professional, sporting or leisure activities), together with the service life of the prosthesis. Patients must also be informed of the need to avoid strenuous effort before attachment and scarring is complete. Repeated and/or excessive stresses may result in the defective attachment of the implant.

Every precaution must be taken to avoid damage to prosthetic components during fitting. Particular attention must be paid to cement particles and/or bone fragments, in the case of cemented prostheses.

In all cases, the prosthesis will not provide the same function as a normal and healthy natural shoulder joint.

Any abnormal pain, unexpected restriction in the amplitude of movement, or variation in the performance of the prosthesis further to attachment and the completion of rehabilitation, must be referred to a medical practitioner.

Components of the MOVE-UP® shoulder prosthesis are not designed for use in conjunction with the components of other manufacturers.

Merge-UP® implants must be fitted using Merge-UP® instruments which are designed for this purpose.

Merge-UP® prosthetic components are single-use components, and must never be reconditioned, resterilized, modified or altered in any way. Any re-implantation of an implant will expose the patient to biological risks which may have severe consequences.

No test for the evaluation of the safety and compatibility of the implant system in a MRI environment (Magnetic Resonance Imaging) has been conducted for the evaluation of associated risks. These include the potential heat-up and/or migration of metallic components. Image artifacts, such as signal drop-out zones and distortions may occur, particularly in the area immediately surrounding the implant.

Any severe incident arising in conjunction with the prosthesis must be reported to Merge-UP and to the competent authority of the country in which the surgeon and/or the patient is established.

## 6 - Associations

Key:

C = Compatible (*diametral mismatch in mm*)

I = Incompatible

NR = Not recommended (*diametral mismatch in mm*)

		Humeral head								
		Ø 37 mm	Ø 39 mm	Ø 41 mm	Ø 43 mm	Ø 45 mm	Ø 47 mm	Ø 49 mm	Ø 51 mm	Ø 53 mm
Cementless humeral short stem	Size 1	C	C	C	NR	NR	NR	NR	NR	NR
	Size 2	NR	NR	NR	C	C	C	NR	NR	NR
	Size 3	NR	NR	NR	NR	NR	NR	C	C	C

		Humeral head, eccentric 1.5 mm								
		Ø 37 mm	Ø 39 mm	Ø 41 mm	Ø 43 mm	Ø 45 mm	Ø 47 mm	Ø 49 mm	Ø 51 mm	Ø 53 mm
ISA™ ONLAY cemented and cementless humeral stems	Size 1	C	C	C	C	C	C	C	C	C
	Size 2	C	C	C	C	C	C	C	C	C
	Size 3	C	C	C	C	C	C	C	C	C
	Size 4	C	C	C	C	C	C	C	C	C
	Size 5	C	C	C	C	C	C	C	C	C
	Size 6	I	C	C	C	C	C	C	C	C
	Size 7	I	I	C	C	C	C	C	C	C
	Size 8	I	I	I	C	C	C	C	C	C
	Size 9	I	I	I	I	C	C	C	C	C

		Humeral head, eccentric 3.5 mm								
		Ø 37 mm	Ø 39 mm	Ø 41 mm	Ø 43 mm	Ø 45 mm	Ø 47 mm	Ø 49 mm	Ø 51 mm	Ø 53 mm
ISA™ ONLAY cemented and cementless humeral stems	Size 1	C	C	C	C	C	C	C	C	C
	Size 2	C	C	C	C	C	C	C	C	C
	Size 3	C	C	C	C	C	C	C	C	C
	Size 4	I	C	C	C	C	C	C	C	C
	Size 5	I	I	C	C	C	C	C	C	C
	Size 6	I	I	I	C	C	C	C	C	C
	Size 7	I	I	I	I	C	C	C	C	C
	Size 8	I	I	I	I	I	C	C	C	C
	Size 9	I	I	I	I	I	I	C	C	C

All ISA™ ONLAY humeral heads, both cemented and cementless, are compatible with all ISA™ ONLAY humeral baseplates.

All ISA™ INLAY humeral stems, both cemented and cementless, are compatible with all ISA™ humeral augments.

All ISA™ INLAY humeral stems, both cemented and cementless, are compatible with the dual taper.

All ISA™ INLAY humeral stems are compatible with all humeral inserts.

All humeral stems for fracture and revision are compatible with all ISA™ augments.

All humeral stems for fracture and revision are compatible with the dual taper.

All humeral stems for fracture and revision are compatible with all distal cortical screws.

All humeral heads are compatible with the dual taper.

		Pegged glenoid/Cemented keeled glenoid			
		Size 1 Reaming Ø 60/80 mm	Size 2 Reaming Ø 60/80 mm	Size 3 Reaming Ø 80/120 mm	Size 4 Reaming Ø 80/120 mm
Humeral head eccentric 1.5 mm and 3.5 mm	Ø 37 mm	C (15.6 mm)	C (20 mm)	NR (24.4 mm)	NR (28.8 mm)
	Ø 39 mm	C (13.6 mm)	C (18 mm)	NR (22.4 mm)	NR (26.8 mm)
	Ø 41 mm	C (11.6 mm)	C (16 mm)	C (20.4 mm)	NR (24.8 mm)
	Ø 43 mm	C (9.6 mm)	C (14 mm)	C (18.4 mm)	NR (22.8 mm)
	Ø 45 mm	NR (7.6 mm)	C (12 mm)	C (16.4 mm)	C (20.8 mm)
	Ø 47 mm	NR (5.6 mm)	C (10 mm)	C (14.4 mm)	C (18.8 mm)
	Ø 49 mm	NR (3.6 mm)	C (8 mm)	C (12.4 mm)	C (16.8 mm)
	Ø 51 mm	NR (1.6 mm)	NR (6 mm)	C (10.4 mm)	C (14.8 mm)
	Ø 53 mm	I	NR (4 mm)	C (8.4 mm)	C (12.8 mm)

		Augmented pegged glenoid			
		8/16° Size 1 right/left	8/16° Size 2 right/left	8/16° Size 3 right/left	8/16° Size 4 right/left
Humeral head eccentric 1.5 mm and 3.5 mm	Ø 37 mm	C (15.6 mm)	C (20 mm)	NR (24.4 mm)	NR (28.8 mm)
	Ø 39 mm	C (13.6 mm)	C (18 mm)	NR (22.4 mm)	NR (26.8 mm)
	Ø 41 mm	C (11.6 mm)	C (16 mm)	C (20.4 mm)	NR (24.8 mm)
	Ø 43 mm	C (9.6 mm)	C (14 mm)	C (18.4 mm)	NR (22.8 mm)
	Ø 45 mm	NR (7.6 mm)	C (12 mm)	C (16.4 mm)	C (20.8 mm)
	Ø 47 mm	NR (5.6 mm)	C (10 mm)	C (14.4 mm)	C (18.8 mm)
	Ø 49 mm	NR (3.6 mm)	C (8 mm)	C (12.4 mm)	C (16.8 mm)
	Ø 51 mm	NR (1.6 mm)	NR (6 mm)	C (10.4 mm)	C (14.8 mm)
	Ø 53 mm	I	NR (4 mm)	C (8.4 mm)	C (12.8 mm)

All ISA™ ONLAY humeral baseplates and all ISA™ humeral augments are compatible with all humeral inserts. The designation of humeral inserts defines the diameters of compatible glenospheres.

All glenospheres are compatible with all glenoid baseplates:

- cementless glenoid baseplates,
- cementless glenoid baseplates with central screw,
- cementless augmented glenoid baseplates with central screw.

All cementless glenoid baseplates are compatible with the peripheral peg and all locked cortical screws and compression cortical screws.

All cementless glenoid baseplates with central screw and all cementless augmented glenoid baseplates with central screw are compatible with the peripheral peg.

All cementless glenoid baseplates with central screw and all cementless augmented glenoid baseplates with central screw are compatible with all locked cortical screws and compression cortical screws.

Central screws are only compatible with glenoid baseplates with a central screw. The selection of the central screw must be consistent with the length of the peg specified in the designations of the latter.

## 7 - Undesirable effects and complications

- Debonding of one or more prosthetic components,
- Displacement of one or more prosthetic components,
- Deformation of polyethylene components,
- Excessive wear of polyethylene components,
- Fracture,
- Infection,
- Pain,
- Luxation or subluxation,
- Contracture,
- Elongation or shortening of the arm,
- Reduction of amplitude of movement.

## 8 - Sterilization

This implant is supplied in a sterile double blister pack, further to sterilization by ionizing gamma radiation with a minimum rating of 25 kGy. The implant must be:

- Kept in its original packaging,
- Stored in a clean and dry location,
- Protected against light, extreme temperatures and sources of ionizing radiation.

It is essential that the expiry date and the integrity of the packaging be checked before opening, in order to detect any trace of deterioration. Damage to the packaging may indicate that the implant is no longer sterile, or has been damaged. In this case, the implant must be returned to MOVE-UP, or destroyed in accordance with regulations in force. The medical device and its packaging are single-use items, and must not be re-used under any circumstances. The medical device and its packaging are not suitable for re-sterilization, either by the user or by a healthcare establishment.

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