


	<i>fr</i>	<i>en</i>	<i>de</i>	<i>es</i>	<i>it</i>	<i>pt</i>
	À usage unique	Single use	Zur einmaligen Verwendung	De un solo uso	Monouso	De utilização única
	Utiliser jusqu'au	Use until	Verwendbar bis	Utilizar hasta el	Utilizzare entro il	A utilizar até
	Référence du catalogue	Catalog reference	Katalogreferenz	Referencia del catálogo	Codice catalogo	Referência do catálogo
	Numéro de lot	Batch code	Chargenbezeichnung	Número de lote	Numero di lotto	Número de lote
	Identifiant unique du dispositif	Unique device identification	Eindeutige Geräte-ID	Identificador único del dispositivo	Identificativo unico del dispositivo	Identificador único do dispositivo
	Dispositif stérilisé par irradiation	Irradiation sterilized device	Sterilisation durch Bestrahlung	Dispositivo esterilizado por irradiación	Dispositivo sterilizzato per irradiazione	Dispositivo esterilizado por irradiação
	Ne pas restériliser	Do not resterilize	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Não esterilizar novamente
	Attention voir notice d'instruction	Caution see IFU	Gebrauchsanweisung beachten	Atención consultar el Manual de instrucciones	Attenzione, consultare il manuale di istruzioni	Atenção, ver o manual de instruções
	Notice d'utilisation	Operator's manual; operating instructions	Gebrauchsanleitung	Manual de utilización	Manuale di istruzioni	Manual de utilização
	Ne pas utiliser si emballage endommagé	Do not use if packaging is damaged	Bei beschädigter Verpackung nicht verwenden	No utilizar si el embalaje está dañado	Non utilizzare se l'imballaggio è danneggiato	Não utilizar se a embalagem estiver avariada
	Crain l'humidité stockage en lieu sec	Dry storage	Vor Nässe schützen, trocken lagern	No tolera la humedad almacenar en un lugar seco	Teme l'umidità, conservare in un luogo asciutto	Teme a humidade, armazenamento em local seco
	Fabricant	Manufacturer	Hersteller	Fabricante	Produttore	Fabricante
	Date de fabrication	Date of manufacture	Herstellungsdatum	Fecha de fabricación	Data di fabbricazione	Data de fabrico
	Fragile	Fragile	Zerbrechlich	Frágil	Fragile	Frágil
	Conserver à l'abri du soleil	Keep away from sunlight	Gegen Sonneinstrahlung geschützt lagern	Conservar protegido del sol	Conservare al riparo dalla luce del sole	Conservar ao abrigo do sol
	Quantité 1	Quantity 1	Menge 1	Cantidad 1	Quantità 1	Quantidade 1
	IRM conditionnel	MR conditional	MRT bedingt	MRI condicional	a compatibilità RM condizionata	IRM condicional
	Dispositif sous ordonnance	Prescription devices	Verschreibungspflichtiges Produkt	Dispositivo dispensado con prescripción médica	Dispositivo su prescrizione medica	Dispositivo sob prescrição médica
	Dispositif médical	Medical device	Medizinprodukt	Dispositivo médico	Dispositivo medico	Dispositivo médico
	Système de barrière stérile double	Double sterile barrier system	System mit doppelter steriler Barriere	Sistema de barrera estéril doble	Sistema de barreira estéril dupla	Doppio sistema di barriera sterile
#	Taille	Size	Größe	Tamaño	Misura	Tamanho
L=	Longueur	Length	Länge	Longitud	Lunghezza	Comprimento
ECC=	Excentré(e)	Eccentric	Exzentrisch	Excéntrico/a	Eccentrico/a	Descentrado(a)
H=	Épaisseur	Thickness	Höhe	Espesor	Spessore	Espessura
RIGHT	Droite	Right	Rechts	Derecho	Destro/a	Direita
LEFT	Gauche	Left	Links	Izquierdo	Sinistro/a	Esquerda
CoCr	Alliage de chrome, cobalt et molybdène	Cobalt Chromium Molybdenum alloy	Chrom-Kobalt-Molybdän-Legierung	Aleación de cromo, cobalto y molibdeno	Lega di cromo, cobalto e molibdeno	Liga de cromo, cobalto e molibdénio
Ti6 Al4 V	Alliage de titane, aluminium et vanadium	Titanium Aluminum Vanadium alloy	Titan-Aluminium-Vanadium-Legierung	Aleación de titanio, aluminio y vanadio	Lega di titanio, alluminio e vanadio	Liga de titânio, alumínio e vanádio
TiN	Revêtement de nitrure de titane	Titanium nitride coating	Titannitrid-Beschichtung	Revestimiento de nitruro de titanio	Rivestimento in nitruro di titanio	Revestimento de nitrato de titânio
Ti	Revêtement de Titane	Titanium coating	Titanbeschichtung	Revestimiento de Titanio	Rivestimento di titanio	Revestimento de Titânio
Ti + HAP	Revêtement de Titane + Hydroxyapatite	Titanium + Hydroxyapatite coating	Titanbeschichtung + Hydroxyapatit	Revestimiento de Titanio + Hidroxiapatita	Rivestimento di titanio + Idrossiapatite	Revestimento de Titânio + Hidroxiapatita

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**1 - Description and intended use**

A partial or total shoulder prosthesis is comprised of various components, all of which are packaged individually and are designed to replace all or part of the natural shoulder joint. A total or partial shoulder arthroplasty is intended to reduce pain and improve the mobility of the patient's joint. The performance capabilities of the device and its service life are closely linked to a multitude of factors, which are dependent upon the patient (physiology, occupation) and upon the surgery undertaken by the surgical team (selection, fitting and fixing of components).

**2 - Materials**

Components	Fixation method		Materials abbreviation	Reference standard
	US	EU		
Ti coated Short Humeral Stem isastemless™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Ti coated Humeral Stem isaonlay™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Uncoated or TIN coated Humeral Stem isaonlay™	CD	CD	Ti6 Al4 V-ELI or Ti6 Al4 V-ELI + TiN	ISO 5832-3
Ti coated Humeral Stem isainlay™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Uncoated or TIN coated Humeral Stem isainlay™	CD	CD	Ti6 Al4 V-ELI or Ti6 Al4 V-ELI + TiN	ISO 5832-3
Ti + HAP coated Humeral Stem isafracture™	CD	CS	Ti6 Al4 V-ELI + Ti + HAP	ISO 5832-3
Ti coated Humeral Stem isarevision™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Dual Taper isafracture™	CD	CS or CD	Ti6 Al4 V-ELI	ISO 5832-3
Uncoated or TIN coated Humeral Baseplate isaonlay™	CD	CS or CD	Ti6 Al4 V-ELI or Ti6 Al4 V-ELI + TiN	ISO 5832-3
Uncoated or TIN coated Humeral Augment isainlay™	CD	CS or CD	Ti6 Al4 V-ELI or Ti6 Al4 V-ELI + TiN	ISO 5832-3
Uncoated or TIN coated Humeral Head isahead™	CD	CS or CD	CoCr or Ti6 Al4 V-ELI + TiN	ISO 5832-12 or ISO 5832-3
Humeral Liner isaliner™	CD	CS or CD	UHMWPE and Ti6 Al4 V-ELI	ISO 5834-1 and ISO 5832-3
Uncoated or TIN coated Glenspheres isasphere™	CD	CS or CD	CoCr or Ti6 Al4 V-ELI + TiN	ISO 5832-12 or ISO 5832-3
Pegged Glenoid isahybrid™	CD	CD	UHMWPE and T40 + Ti	ISO 5834-1 and ISO 5832-2
Cemented Keeled Glenoid isalegacy™	CD	CD	UHMWPE and T40	ISO 5834-1 and ISO 5832-2
Augmented Pegged Glenoid isahybrid+™	CD	CD	UHMWPE and T40 + Ti	ISO 5834-1 and ISO 5832-2
Ti coated Humeral baseplate isareverse™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Ti coated Glenoid Baseplate with Central Screw isareverse™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Ti coated Augmented Glenoid Baseplate with Central Screw isareverse+™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Ti coated Locked Peripheral Peg isareverse™	CD	CS	Ti6 Al4 V-ELI + Ti	ISO 5832-3
Cortical Screws, Central Screws, Obturator	CD	CS or CD	Ti6 Al4 V-ELI	ISO 5832-3

CD = Cemented; CS = Cementless

**3 - Indications**

	Intact glenoid humeral hemiarthroplasties	Total shoulder arthroplasties performed with anatomic prostheses	Total shoulder arthroplasties performed with inverted prostheses
→ Components ↓	Humeral head, screws	Glenoid, Humeral head, Dual taper	Glenoid baseplate, Glensphere, Humeral liner, Humeral baseplate, Humeral Augment, screws
Stemless stems	- Primary osteoarthritis - Post-traumatic osteoarthritis - Aseptic osteonecrosis	- Gleno humeral osteoarthritis - Primary gleno humeral osteoarthritis - Post-traumatic osteoarthritis	
Onlay stems		- Gleno humeral osteoarthritis - Rheumatoid arthritis - Osteoarthritis with glenoid posterior wear - Osteoarthritis with B2 glenoid	- Primary Gleno humeral osteoarthritis - Gleno humeral osteoarthritis grade 1 to 5 - Cuff tear arthropathy - Massive rotator cuff tear - Post-traumatic osteoarthritis - Avascular necrosis - Chronic instability - Glenoid erosion E0 to E3 - Rheumatoid arthritis
Inlay stems		- Primary gleno humeral osteoarthritis	- Gleno humeral osteoarthritis - Primary rotator cuff deficiency - Failure of rotator cuff operation - Failure of arthroplasty - Prosthetic instability - Proximal humerus malunion
Fracture stems	- Humeral head displacement - Three-part fractures - Four-part fractures - Proximal humeral fracture-dislocations - Comminuted fractures of the humeral neck		- Acute proximal humeral fractures - Fractures of the humeral neck - Glenoid defects
Revision stems			- Gleno humeral osteoarthritis - Cuff tear arthropathy - Severe rotator cuff deficiency - Failure of rotator cuff operation - Failure of hemiarthroplasty - Failure of total arthroplasty - Humeral fracture - Proximal humeral fracture and sequelae - Shoulder pseudoparalysis - Proximal humerus malunion

Titanium Nitride coated Ti6 Al4 V-ELI (Ti6 Al4 V-ELI + TiN) components are intended for allergic patients or for patients with suspected sensitivity to CoCr alloy constituents: nickel, manganese, chromium, cobalt and/or molybdenum.

#### 4 - Counter-indications

The implantation of a total or partial shoulder prosthesis, in common with any surgical procedure, has attendant risks, specifically associated with anaesthesia and with the operation itself.

##### General risks:

- Systemic or local infections,
- Severe vascular or neurological deficiencies
- Bone loss or poor bone quality, which is liable to compromise the stability of the implant
- Allergies to materials present in the implant
- Bone tumour in the region of the implant

##### Anatomical prostheses:

- Severe deficiency of the cuff
- Insufficient bone stock for the attachment of implants
- Permanent instability

#### 5 - Warnings and precautions

Prior to any clinical use, the surgeon must be fully versed in all aspects of the operation, the implants and associated instruments together with their applications and limits.

Surgeons must notify their patients of the limitations of the prosthesis, and of the associated impact upon their activities (professional, sporting or leisure activities), together with the service life of the prosthesis. Patients must also be informed of the need to avoid strenuous effort before attachment and scarring is complete.

Repeated and/or excessive stresses may result in the defective attachment of the implant.

Every precaution must be taken to avoid damage to prosthetic components during fitting. Particular attention must be paid to cement particles and/or bone fragments.

In all cases, the prosthesis will not provide the same function as a normal and healthy natural shoulder joint.

Any abnormal pain, unexpected restriction in the amplitude of movement, or variation in the performance of the prosthesis further to attachment and the completion of rehabilitation, must be referred to a medical practitioner.

Each humeral stem box contains an implant card to complete by sticking the adhesive labels of the implanted components. This card contains patient information.

Components of the isa™ MOVE-UP shoulder prosthesis are not designed for use in conjunction with the components of other manufacturers.

MOVE-UP implants must be fitted using MOVE-UP instruments which are designed for this purpose.

MOVE-UP prosthetic components are single-use components, and must never be reconditioned, resterilized, modified or altered in any way. Any re-implantation of an implant will expose the patient to biological risks which may have severe consequences.

Non-clinical testing has demonstrated that the MOVE-UP devices are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 or 3.0 Tesla,
- Horizontal bore MRI system,
- Maximum spatial gradient magnetic field of 1900 gauss/cm (19T/m),
- $B0 \cdot |dB0/dr|$  product lower or equal to 48T<sup>2</sup>/m,
- RF whole body transmit/receive coil use only,

- Whole body averaged SAR (Specific Absorption Rate) limited to Normal operating mode (WB-SAR  $\leq 2$ W/kg).

MRI related torque and displacement force: No restriction for the position with 1.5T and 3.0T.

MRI related heating:

	1.5 Tesla (64MHz)	3.0 Tesla (128MHz)
Maximum MR system, whole body averaged SAR (WB-SAR)	2 W/kg after 15 minutes of continuous scanning	2 W/kg after 15 minutes of continuous scanning
Calorimetry measured values, whole body averaged SAR (WB-SAR)	2.08 +/- 0.81 W/kg	2.11 +/- 0.88W/kg
Maximal temperature increase measured on the device	9.8 +/- 1.0°C	4.9 +/- 1.0°C

MRI related artifact:

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the medical device. MOVE-UP has not evaluated the size of artifact.

Any severe incident arising in conjunction with the prosthesis must be reported to MOVE-UP and to the competent authority of the country in which the surgeon and/or the patient is established.

#### 6 - Associations

See table of associations.

Key:

C = Compatible (*diametral mismatch in mm*)

I = Incompatible

NR = Not recommended (*diametral mismatch in mm*)

All onlay humeral stems, both cemented and cementless, are compatible with all onlay humeral baseplates.

All inlay humeral stems, both cemented and cementless, are compatible with all humeral augments.

All inlay humeral stems, both cemented and cementless, are compatible with the dual taper.

All inlay humeral stems are compatible with all humeral inserts.

All humeral stems for fracture and revision are compatible with all augments.

All humeral stems for fracture and revision are compatible with the dual taper.

All humeral stems for fracture and revision are compatible with all distal cortical screws.

All humeral stems for fracture and revision are compatible with all humeral inserts.

All humeral heads are compatible with the dual taper.

All onlay humeral baseplates and all humeral augments are compatible with all humeral inserts.

The designation of humeral inserts defines the diameters of compatible glenospheres.

All glenospheres are compatible with all glenoid baseplates:

- cementless glenoid baseplates,
- cementless glenoid baseplates with central screw hole,
- cementless augmented glenoid baseplates with central screw hole.

All glenoid baseplates are compatible with the peripheral peg and all locked cortical screws and compression cortical screws.

All glenoid baseplates with central screw hole are compatible with all central screws and the obturator.

**7 - Undesirable effects and complications**

- Debonding of one or more prosthetic components
- Displacement of one or more prosthetic components
- Disassembly of one or more prosthetic components
- Rupture of one or more prosthetic components
- Deformation of prosthetic components
- Excessive wear of prosthetic components
- Fracture
- Infection
- Pain
- Luxation or subluxation
- Contracture
- Elongation or shortening of the arm
- Reduction of movement amplitude (impingement)
- Bone remodeling (notching)
- Allergies to materials present in the implant

**8 - Sterilization**

This implant is supplied in a sterile double blister pack, further to sterilization by ionizing gamma radiation with a minimum rating of 25 kGy. The implant must be:

- Kept in its original packaging,
- Stored in a clean and dry location,
- Protected against light, extreme temperatures and sources of ionizing radiation.

It is essential that the expiry date and the integrity of the packaging be checked before opening, in order to detect any trace of deterioration. A damaged or inadvertently opened packaging is not fit for use. In this case, the implant must be returned to MOVE-UP, or destroyed in accordance with regulations in force. The medical device and its packaging are single-use items and must not be re-used under any circumstances. The medical device and its packaging are not suitable for re-sterilization, either by the user or by a healthcare establishment.

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